












钛特利® T-I-T-L-I-S Implant System TL Product List Surgery		Tissue Level					
		Φ 4.1 mm		Φ 4.7 mm			
							
Dental Implant	SLA	020.4148.06	021.4148.06	020.4748.06	021.4748.06	020.4760.06	021.4760.06
	6.5 mm	020.4148.08	021.4148.08	020.4748.08	021.4748.08	020.4760.08	021.4760.08
	8 mm	020.4148.10	021.4148.10	020.4748.10	021.4748.10	020.4760.10	021.4760.10
	10 mm	020.4148.12	021.4148.12	020.4748.12	021.4748.12	020.4760.12	021.4760.12
	12 mm	020.4148.14	021.4148.14	-	-	-	-
Dental Drill (3-blade Drill Sharp/Round Drill)	Φ 1.8 mm	Φ 1.5 mm		Φ 2.3 mm		Φ 3.5 mm	
							
Dental Drill (Pilot Drill)		180.1800.33		181.1500.33		181.2300.33	
		181.1500.33		181.2300.33		181.3500.33	
Dental Drill (Final Drill)	Φ 2.0 mm	182.2000.33		182.2000.41			
		182.2000.33		182.2000.41			
Dental Drill (Final Drill)	Φ 2.8 mm	Φ 3.4 mm		Φ 4.0 mm			
							
Alignment Pin /Depth Gauge		183.2800.33		183.3400.33		183.4000.33	
		183.2800.41		183.3400.41		183.4000.41	
Dental Drill (Profile Drill)	Φ 1.9 mm	Φ 2.7 mm		Φ 3.3 mm		Φ 3.9 mm	
							
Dental Drill (Tap Drill)		190.1900.28		191.2700.26		191.3300.26	
		191.2700.26		191.3300.26		191.3900.26	
Dental Drill (Tap Drill)	Φ 4.0 mm	Φ 4.6 mm					
							
Dental Drill (Tap Drill)		184.4000.26		184.4000.33		184.4600.26	
		184.4000.33		184.4600.26		184.4600.33	
Dental Drill (Tap Drill)	Φ 4.1 mm	Φ 4.7 mm					
							
Dental Drill (Tap Drill)		185.4100.23		185.4100.33		185.4100.41	
		185.4100.33		185.4100.41		185.4700.23	

钛特利® T-I-T-L-I-S Implant System TL Product List Surgery		Tissue Level					
		Φ 4.1 mm		Φ 4.7 mm			
							
Dental Implant	SLA	020.4148.06	021.4148.06	020.4748.06	021.4748.06	020.4760.06	021.4760.06
	6.5 mm	020.4148.08	021.4148.08	020.4748.08	021.4748.08	020.4760.08	021.4760.08
	8 mm	020.4148.10	021.4148.10	020.4748.10	021.4748.10	020.4760.10	021.4760.10
	10 mm	020.4148.12	021.4148.12	020.4748.12	021.4748.12	020.4760.12	021.4760.12
	12 mm	020.4148.14	021.4148.14	-	-	-	-
Closure Screw	Φ 3.4 mm	Φ 4.4 mm					
							
Healing Cap	H 0.5 mm	H 0.5 mm					
	170.0001.00	170.0002.00					
Healing Cap	Φ 5.3 mm	Φ 6.7 mm					
							
Healing Cap	H 1 mm	H 2.5 mm	H 4 mm	H 1 mm	H 2.5 mm	H 4 mm	
	040.5300.10	040.5300.25	040.5300.40	040.6700.10	040.6700.25	040.6700.40	

Regular Neck Wide Neck

Sino - Swiss Joint Production

Quality and customer
satisfaction are our
eternal pursuit










Wechat Public
Account








Rui Qi Medical Technology (Jiaxing) Co., Ltd.



Titlis Implant System Product List

钛特利® T-I-T-L-I-S Implant System TL Product List Prosthetics		Tissue Level					
		Φ 4.1 mm		Φ 4.7 mm			
							
Dental Implant	SLA	020.4148.06	021.4148.06	020.4748.06	021.4748.06	020.4760.06	021.4760.06
	6.5 mm	020.4148.08	021.4148.08	020.4748.08	021.4748.08	020.4760.08	021.4760.08
	8 mm	020.4148.10	021.4148.10	020.4748.10	021.4748.10	020.4760.10	021.4760.10
	10 mm	020.4148.12	021.4148.12	020.4748.12	021.4748.12	020.4760.12	021.4760.12
	12 mm	020.4148.14	021.4148.14	-	-	-	-
Abutment	Straight Abutment (Screw Retained)	Straight Abutment (Screw Retained)					
	Φ 3.4 mm	Φ 4.4 mm					
Abutment							
	AH 2 mm	AH 2 mm					
Abutment	Straight Abutment (Cement Retained)	Straight Abutment (Cement Retained)					
	Φ 3.4 mm	Φ 4.4 mm					
Abutment							
	AH 5.5 mm	AH 5.5 mm					
Abutment	061.3400.55	061.4400.55					
	Angled Abutment, 20° (Cement Or Screw Retained)	Angled Abutment, 20° (Cement Or Screw Retained)					
Abutment	Φ 3.4 mm	Φ 4.4 mm					
							
	AH 5.5 mm	AH 5.5 mm					
Abutment	062.3400.55	062.4400.55					
	Φ 5.2 mm	Φ 6.5 mm					
Impression Post							
	L 11 mm	L 15 mm	L 11 mm	L 10.7 mm	L 14.7 mm	L 10.7 mm	
Impression Post	120.5200.13	122.5200.13	121.5200.15	120.6500.13	122.6500.13	121.6500.15	
Implant Analog	Φ 11.8 mm	Φ 13.8 mm					
	080.4148.12	080.4760.14					
Plastic Coping							
	H 10 mm	H 10 mm	H 10 mm	H 10 mm	H 10 mm	H 10 mm	
Plastic Coping	303.5100.10	304.5100.10	305.5100.10	306.5100.10	303.6600.10	304.6600.10	305.6600.10
Occlusal Screw							
	L 5.2 mm	L 4 mm					
Occlusal Screw	220.2100.52	221.3200.40					

钛特利® T-I-T-L-I-S Implant System Digitalized Product List Prosthetics		Tissue Level					
		Φ 4.1 mm		Φ 4.7 mm			
							
Dental Implant	SLA	020.4148.06	021.4148.06	020.4748.06	021.4748.06	020.4760.06	021.4760.06
	6.5 mm	020.4148.08	021.4148.08	020.4748.08	021.4748.08	020.4760.08	021.4760.08
	8 mm	020.4148.10	021.4148.10	020.4748.10	021.4748.10	020.4760.10	021.4760.10
	10 mm	020.4148.12	021.4148.12	020.4748.12	021.4748.12	020.4760.12	021.4760.12
	12 mm	020.4148.14	021.4148.14	-	-	-	-
Ti Base							
		For Crown			For Crown		
Scanbody		HPLA-0311 (GB-B_TTL-TL-48)			HPLA-0311 (GB-B_TTL-TL-60)		
							
Digital Implant Analog		SC-B-4.0 (SC-B-4.0_TTL-TL-48)			SC-B-4.5 (SC-B-4.5_TTL-TL-60)		
							
CAD/CAM Abutment		AN-D/4.0 (AN-D/4.0_TTL-TL-48)			AN-D/4.5 (AN-D/4.5_TTL-TL-60)		
							

 钛特利® T-I-T-L-I-S Implant System BL Product List Surgery		Bone Level								
		Φ 3.5 mm			Φ 4.1 mm			Φ 4.7 mm		
										
		Φ 3.5 mm			Φ 4.1 mm			Φ 4.7 mm		
Dental Implant	SLA									
	8 mm	010.3500.08			010.4100.08			010.4700.08		
	10 mm	010.3500.10			010.4100.10			010.4700.10		
	12 mm	010.3500.12			010.4100.12			010.4700.12		
	14 mm	010.3500.14			010.4100.14			010.4700.14		
Dental Drill (3-blade Drill Sharp/Round Drill)	Φ 1.8 mm		Φ 1.5mm		Φ 2.3 mm		Φ 3.5 mm			
										
		180.1800.33		181.1500.33		181.2300.33		181.3500.33		
Dental Drill (Pilot Drill)	Φ 2.0 mm									
										
		182.2000.33				182.2000.41				
Dental Drill (Final Drill)	Φ 2.8 mm		Φ 3.4 mm		Φ 4.0 mm					
										
	183.2800.33		183.3400.33		183.4000.33					
		183.2800.41		183.3400.41		183.4000.41				
Alignment Pin /Depth Gauge	Φ 1.9 mm		Φ 2.7 mm		Φ 3.3 mm		Φ 3.9 mm			
										
		190.1900.28		191.2700.26		191.3300.26		191.3900.26		
Dental Drill (Profile Drill)	Φ 3.4 mm		Φ 4.0 mm		Φ 4.6 mm					
										
		184.3400.26		184.3400.33		184.4000.26		184.4600.33		
Dental Drill (Tap Drill)	Φ 3.5 mm			Φ 4.1 mm			Φ 4.7 mm			
										
		185.3500.23	185.3500.33	185.3500.41	185.4100.23	185.4100.33	185.4100.41	185.4700.23	185.4700.33	185.4700.41